


Mobile-App to evaluate tumor visualization strategies in Augmented Reality

by Felix Meyer and Remke Albrecht






NATIONAL CENTER
FOR TUMOR DISEASES
PARTNER SITE DRESDEN
UNIVERSITY CANCER CENTER UCC


Supported by:
German Cancer Research Center
University Hospital Carl Gustav Carus Dresden
Carl Gustav Carus Faculty of Medicine, TU Dresden
Helmholtz-Zentrum Dresden-Rossendorf

Environment

- 3D printed **silicon kidney** of the phantom
- 3D models of the **kidney, calyx** and **tumor**
- The app should run on portable devices




Tasks

- 
- Task 1 - Tracking
 - Task 2 - Visualization
 - Task 3 - Graphical user interface
 - Task 4 - Track user input and used settings
 - Evaluation of usability


Task 1 - Tracking

- Vuforia Augmented Reality SDK
 - Kidney with simple geometry (few possible feature points)
 - No texture information
 - 3D printing artefacts
- Challenges:
 - Kidney with simple geometry (few possible feature points)
 - No texture information
 - 3D printing artefacts
- Solutions:
 - Use multiple hand picked target views for the object
 - Not too many to minimize performance overhead
 - Maximize visibility with external light sources or integrated flashlight
 - Interpolate tracking with acceleration sensor



Task 2 - Visualization

- **Unity** as engine:
 - Integrated **Vuforia** support
 - **Simple** and **adjustable** render pipeline
 - **Android** deployment
 - Portable performance
- Big **new problem** with transparency:
 - Only tested with single transparent and multiply opaque objects
 - Multiple transparent objects require new rendering technique






Task 2 - Visualization

The Problem:

- Happens between **transparent objects**
- Happens with single **concave** transparent object
- Unitys solution for multiple transparent objects needs a **clearly defined object order**
 - Solves only transparency problem between objects and only if they **don't intersect**
 - We have **concave objects** and **intersecting objects**

→ Won't work for us



Task 2 - Visualization


Tested transparency solutions:

Technique	Pro	Con
Order dependent transparency:		
Object sorting	<ul style="list-style-type: none">• Solves some transparency issues between objects	<ul style="list-style-type: none">• Won't work with intersecting objects• Concave objects still a problem
Vertex sorting	<ul style="list-style-type: none">• Solves all transparency issues	<ul style="list-style-type: none">• Difficult implementation• Far too heavy for mobile devices
Order independent transparency:		
Weighted transparency	<ul style="list-style-type: none">• Works well with very transparent objects• Average blend between transparent surfaces• Fast	<ul style="list-style-type: none">• intersections between surfaces are invisible• falls apart for opaque objects
Depth peeling	<ul style="list-style-type: none">• Trade-off quality and rendertime• Solves all transparency issues• works good on mobile devices	<ul style="list-style-type: none">• Maximum number of overlaying transparent layers

Task 2 - Visualization


Final solution: **Depth peeling**

- Gives best control over performance and quality
- Object intersections are visible even on transparent objects
- Currently use 6 transparency layers
 - after 6 overlapping transparent vertices, artefacts are theoretically visible (not visible in practice)
 - gives acceptable performance on tested tablet and phones (Renders scene 6 times)




Task 3 - Graphical User Interface

- Previous solution
 - Select Object
 - Select visual channel
 - Select value the visual channel should take
- New solution
 - Enter name
 - Select one of six presets
 - Reselect until user is satisfied
 - Fine tune visual channels
 - Submit




Step 1: Enter name



Step 2: Select Preset



Step 3: Confirm




Step 4: Select Object

Step 5: Select Visual Channel

Step 6: Select Value

Step 7: Submit




Back to Presets

Preset Comparison

Redo

Undo





Task 4 - Track user input and used settings

- Save **user name**
- **Track** every action with timestamp
- Save last set settings

E.g. History_User1.csv

Object	Effect	Value	Time
Niere	_Alpha	0,6091537	23,88399
Niere	_Glossiness	6,119171	27,31732
Tumor	_Glossiness	9,075993	32,00634
Tumor	_Alpha	0,819171	36,37442



Kidney

Transparency
Falloff

Back To Presets




vuforia™

Undo

Submit

Redo



Thank you for your attention!